

**Nucleic acid adduct for in situ hybridisation to give nucleic acid sequences - comprises water soluble basic polymer and nucleic acid sequence formed by polymerase chain reaction and/or reverse transcriptase reaction**

**Patent number:** DE4038293  
**Publication date:** 1992-06-04  
**Inventor:** SCHLECHTE HORST DR (DE); WESTERMANN PETER DR (DE)  
**Applicant:** INST MOLEKULARBIOLOGIE AK (DE); ZENTRALINSTITUT FUER KREBSFORS (DE)  
**Classification:**  
 - international: C07H21/04; C12P19/34  
 - european: C12N15/10, C12Q1/68B14  
**Application number:** DE19904038293 19901128  
**Priority number(s):** DE19904038293 19901128

**Abstract of DE4038293**

Nucleic acid adduct comprises (a) a water soluble, basic polymer, adsorbed on cells in histological or cytological preps.; and (b) a nucleic acid sequence contained in the cells and is formed by a polymerase chain reaction and/or reverse transcriptase reaction. Adducts are amplifiable and there is a bond between the two components. Basic polymer is pref. polyethylene-imine or DEAE-dextran. **USE/ADVANTAGE** - Can be used for in situ hybridisation to give special nucleic acid sequences. Double-stranded or single-stranded DNA or RNA, of viral or cellular origin, is bound in situ to the basic polymer and the adduct is amplified by a polymerase chain reaction and hybridised with a labelled gene probe. In an example, cervix carcinoma tissue was fixed with neutral formaldehyde soln. and embedded in paraffin. A slice of this was placed on a microscopic slide which had been pre-treated with 3-(triethoxysilyl)-propylamine. This was incubated with 5% polyethylene imine (mol.wt. 30,000-40,000) in water for 30 mins. at 37 deg.C. Glass carrier was rinsed with water and dried. Cover slide was placed over the tissue and a conventional polymerase chain reaction was carried out by dipping the prepn. in a thermostatic bath, followed by hybridisation with biotin labelled human papilloma virus type 16-DNA. (Dwg.0/0)

---

Data supplied from the **esp@cenet** database - Worldwide

**Caren Burgoon**

---

**From:** Simon.Kiddle@Mewburn.com  
**Sent:** Monday, February 23, 2004 10:08 AM  
**To:** Pat Hagan  
**Subject:** Document request for DE 4038293/IDS for Matthew Baker Charge Switch applications

Hi Pat

This is to follow up on the phone call from your assistant Karen.

You can find an abstract for De 4038293 at:

<http://v3.espacenet.com/textdoc?DB=EPODOC&IDX=DE4038293>

I believe that this should be sufficient for the US PTO, but let me know if you need anything further.

Best regards.

Simon

**MEWBURN ELLIS**

Phone: +44 (0)117 926 6411

Fax: +44 (0)117 926 5692

London Fax: +44 (0)20-7240-9339

Website: <http://www.mewburn.com>

The information in this e-mail message may be privileged and is confidential information intended only for the use of the recipient named above. If the reader of this message is not the intended recipient, or the employee or agent responsible to deliver it to the intended recipient, any use, distribution or copying of this communication is strictly prohibited. If you have received this communication in error, please notify us immediately by telephone and destroy the original message from your electronic files.

02/23/2004